



Effective Date:  
Monday, December 02, 2019

## Test Updates

In our continuing effort to provide you with the highest quality toxicology laboratory services available, we have compiled important changes regarding a number of tests we perform. Listed below are the types of changes that may be included in this notification, effective Monday, December 02, 2019

---

**Test Changes** - Tests that have had changes to the method/ CPT code, units of measurement, scope of analysis, reference comments, or specimen requirements.

**Discontinued Tests** - Tests being discontinued with alternate testing suggestions.

---

Please use this information to update your computer systems/records. These changes are important to ensure standardization of our mutual laboratory databases.

If you have any questions about the information contained in this notification, please call our Client Support Department at (866) 522-2206. Thank you for your continued support of NMS Labs and your assistance in implementing these changes.

The CPT Codes provided in this document are based on AMA guidelines and are for informational purposes only. NMS Labs does not assume responsibility for billing errors due to reliance on the CPT Codes listed in this document.



Effective Date:  
Monday, December 02, 2019

## Test Updates

Test Code	Test Name	Test Name	Method / CPT Code	Specimen Req.	Stability	Scope	Units	Reference Comments	Discontinue
8890OF	Amphetamines Panel (Qualitative), Oral Fluid (Saliva)		•						
8889OF	DUID-DRE Panel (Qualitative) ProofPOSITIVE®, Oral Fluids (CSA)		•						
8898OF	Drugs of Abuse (6 Panel) (Qualitative), Oral Fluid (Saliva)		•						
8897OF	Drugs of Abuse (7 Panel) (Qualitative), Oral Fluid (Saliva)		•						
9176U	Fentanyl Screen, Urine				•				
10037U	Fentanyl and Metabolite (Qualitative) (DUID/DRE), Urine (Forensic) (CSA)				•				
54142U	Fentanyl and Metabolite Confirmation (Qualitative) (DUID/DRE), Urine				•				
5640U	Fentanyl and Metabolite Confirmation, Urine				•				
55041U	Fentanyl and Metabolite Confirmation, Urine (CSA)				•				
2079U	Fentanyl and Metabolite, Urine				•				



## Test Updates

### Test Changes

#### 8890OF Amphetamines Panel (Qualitative), Oral Fluid (Saliva)

Summary of Changes: Methods/CPT Codes were changed [LC-MS/MS (80324, 80359)]

Scope of Analysis: LC-MS/MS (80324, 80359): Amphetamine, Methamphetamine, MDA, MDMA  
Method (CPT Code)

#### 8889OF DUID-DRE Panel (Qualitative) ProofPOSITIVE®, Oral Fluids (CSA)

Summary of Changes: Methods/CPT Codes were changed [LC-MS/MS (80324,80346,80353,80358,80359, 80361, 80362, 80365, 83992)]

Scope of Analysis: LC-MS/MS (80324,80346,80353,80358,80359, 80361, 80362, 80365, 83992):  
Method (CPT Code) Amphetamine, Methamphetamine, MDA, MDMA, Diazepam, Nordiazepam, Oxazepam, Temazepam, Chlordiazepoxide, Lorazepam, Clonazepam, Alprazolam, Midazolam, Codeine - Free, Morphine - Free, Hydrocodone - Free, 6-MAM - Free, Hydromorphone - Free, Oxycodone - Free, Oxymorphone - Free, Dihydrocodeine - Free, Cocaine, Benzoyllecgonine, Cocaethylene, Methadone, EDDP, Phencyclidine, Dextromethorphan

#### 8898OF Drugs of Abuse (6 Panel) (Qualitative), Oral Fluid (Saliva)

Summary of Changes: Methods/CPT Codes were changed [LC-MS/MS (80324,80346,80353,80358,80359, 80361, 80362, 80365, 83992)]

Scope of Analysis: LC-MS/MS (80324,80346,80353,80358,80359, 80361, 80362, 80365, 83992):  
Method (CPT Code) Amphetamine, Methamphetamine, MDA, MDMA, Diazepam, Nordiazepam, Oxazepam, Temazepam, Chlordiazepoxide, Lorazepam, Clonazepam, Alprazolam, Midazolam, Cocaine, Benzoyllecgonine, Cocaethylene, Methadone, EDDP, Codeine - Free, Morphine - Free, Hydrocodone - Free, 6-MAM - Free, Hydromorphone - Free, Oxycodone - Free, Oxymorphone - Free, Dihydrocodeine - Free, Phencyclidine, Dextromethorphan

#### 8897OF Drugs of Abuse (7 Panel) (Qualitative), Oral Fluid (Saliva)

Summary of Changes: Methods/CPT Codes were changed [LC-MS/MS (80324,80346,80353,80358,80359, 80361, 80362, 80365, 83992)]

Scope of Analysis: LC-MS/MS (80324,80346,80353,80358,80359, 80361, 80362, 80365, 83992):  
Method (CPT Code) Amphetamine, Methamphetamine, MDA, MDMA, Diazepam, Nordiazepam, Oxazepam, Temazepam, Chlordiazepoxide, Lorazepam, Clonazepam, Alprazolam, Midazolam, Cocaine, Benzoyllecgonine, Cocaethylene, Methadone, EDDP, Codeine - Free, Morphine - Free, Hydrocodone - Free, 6-MAM - Free, Hydromorphone - Free, Oxycodone - Free, Oxymorphone - Free, Dihydrocodeine - Free, Phencyclidine, Dextromethorphan  
LC-MS/MS (80349): Delta-9 THC

#### 9176U Fentanyl Screen, Urine

Summary of Changes: Stability was changed.



Effective Date:  
Monday, December 02, 2019

## Test Updates

### Test Changes

Stability: Room Temperature: 1 month(s)  
Refrigerated: 1 month(s)  
Frozen (-20 °C): 15 month(s)

#### **10037U Fentanyl and Metabolite (Qualitative) (DUID/DRE), Urine (Forensic) (CSA)**

Summary of Changes: Stability was changed.

---

Stability: Room Temperature: 1 month(s)  
Refrigerated: 1 month(s)  
Frozen (-20 °C): 15 month(s)

#### **54142U Fentanyl and Metabolite Confirmation (Qualitative) (DUID/DRE), Urine**

Summary of Changes: Stability was changed.

---

Stability: Room Temperature: 1 month(s)  
Refrigerated: 1 month(s)  
Frozen (-20 °C): 15 month(s)

#### **55041U Fentanyl and Metabolite Confirmation, Urine (CSA)**

Summary of Changes: Stability was changed.

---

Stability: Room Temperature: 1 month(s)  
Refrigerated: 1 month(s)  
Frozen (-20 °C): 15 month(s)

#### **5640U Fentanyl and Metabolite Confirmation, Urine**

Summary of Changes: Stability was changed.

---

Stability: Room Temperature: 1 month(s)  
Refrigerated: 1 month(s)  
Frozen (-20 °C): 15 month(s)

#### **2079U Fentanyl and Metabolite, Urine**

Summary of Changes: Stability was changed.

---

Stability: Room Temperature: 1 month(s)  
Refrigerated: 1 month(s)  
Frozen (-20 °C): 15 month(s)