



NMS Labs

CONFIDENTIAL

200 Welsh Road, Horsham, PA 19044-2208
Phone: (215) 657-4900 Fax: (215) 657-2972
e-mail: nms@nmslabs.com

Robert A. Middleberg, PhD, F-ABFT, DABCC-TC, Laboratory Director

Demo Report

Report Issued 10/18/2018 11:30
Last Report Issued 11/02/2010 08:26

88888
Clinical Example Report
Attn: IT Department
200 Welsh Road
Horsham, PA 19044-2208

Patient Name 1272SP-POS
Patient ID 1272SP-POS
Chain 10002233
Age Not Given **DOB** Not Given
Gender Not Given
Workorder 10002233
Received 11/02/2010 08:21

Sample ID 10002233-001
Matrix Serum or Plasma
Patient Name 1272SP-POS
Patient ID 1272SP-POS
Container Type Clear vial

Collect Dt/Tm Not Given
Source Not Given

Approx Vol/Weight Not Given

Receipt Notes None Entered

Analysis and Comments	Result	Units	Reporting Limit	Notes
-----------------------	--------	-------	-----------------	-------

1272SP Citalopram, Serum/Plasma

Analysis by High Performance Liquid Chromatography/
Tandem Mass Spectrometry (LC-MS/MS)

Citalopram / Escitalopram	50	ng/mL	5.0	
Synonym(s): Celexa®/Lexapro®				

Steady-state serum or plasma levels from patients on a daily regimen of 30 to 60 mg Citalopram: 9 - 200 ng/mL.

This test is not Chiral specific. Patients who have taken Escitalopram (Lexapro®), as opposed to Racemic Citalopram (Celexa®), within the past 3 days may have falsely elevated values.