



**NMS Labs**

**CONFIDENTIAL**

200 Welsh Road, Horsham, PA 19044-2208  
Phone: (215) 657-4900 Fax: (215) 657-2972  
e-mail: nms@nmslabs.com

Robert A Middleberg, PhD, F-ABFT, DABCC-TC, Laboratory Director

**Demo Report**

**Report Issued** 09/09/2011 14:46

**Last Report Issued** 11/02/2010 08:27

**88888**

Clinical Dummy reports  
Attn: IT Department  
200 Welsh Road  
Horsham, PA 19044-2208

**Patient Name** 1965U-POS

**Patient ID** 1965U-POS

**Chain** 10002227

**Age** Not Given **Gender** Not Given

**Workorder** 10002227

**Received** 11/02/2010 07:19

**Sample ID** 10002227-001  
**Matrix** Urine  
**Patient Name** 1965U-POS  
**Patient ID** 1965U-POS  
**Container Type** Clear vial

**Collect Dt/Tm** Not Given  
**Source** Not Given

**Approx Vol/Weight** Not Given

**Receipt Notes** None Entered

Analysis and Comments	Result	Units	Reporting Limit	Notes
-----------------------	--------	-------	-----------------	-------

**1965U Escitalopram, Urine**

Analysis by High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)

Citalopram / Escitalopram	50	ng/mL	5.0	
Synonym(s): Celexa®/Lexapro®				

This test is not Chiral specific. Patients who have taken Racemic Citalopram (Celexa®), as opposed to Escitalopram (Lexapro®), within the past 3 days may have falsely elevated values.