



NMS Labs

CONFIDENTIAL

200 Welsh Road, Horsham, PA 19044-2208
Phone: (215) 657-4900 Fax: (215) 657-2972
e-mail: nms@nmslabs.com

Robert A. Middleberg, PhD, F-ABFT, DABCC-TC, Laboratory Director

Demo Report

Report Issued 04/21/2022 17:40
Last Report Issued 01/28/2013 09:41

88888
Clinical Example Report
Attn: Example Reports
200 Welsh Road
Horsham, PA 19044

Patient Name 2492U
Patient ID 2492U
Chain 13000442
DOB Not Given
Sex Not Given
Workorder 13000442
Received 01/28/2013

Lab ID 13000442-001
Matrix Urine
Patient Name 2492U
Patient ID 2492U
Container Type Clear vial

Collect Dt/Tm Not Given
Source Not Given

Approx Vol/Weight Not Given

Receipt Notes None Entered

Analysis and Comments	Result	Units	Reporting Limit	Notes
2492U Lead, Urine				
Analysis by Colorimetry (C)				
Creatinine	1000	mg/L	100	
U.S. Population (10th - 90th percentiles, median) All participants: 335 - 2370 mg/L, median 1180 (n=22,245) Males: 495 - 2540 mg/L, median 1370 (n=10,610) Females: 273 - 2170 mg/L, median 994 (n=11,635)				
Analysis by Inductively Coupled Plasma/Mass Spectrometry (ICP/MS)				
Lead	None Detected	mcg/L	5.0	
Normally: less than 5 mcg/L. Various states require that levels above certain cutoffs must be reported to the state in which the patient resides. Please contact NMS Labs if you need assistance in supplying your state with the required information.				
Lead (Creatinine corrected)	None Detected	mcg/g Creat	5.0	

Results for sample 13000442-001 are continued on next page



NMS Labs

CONFIDENTIAL

200 Welsh Road, Horsham, PA 19044-2208

Phone: (215) 657-4900 Fax: (215) 657-2972

e-mail: nms@nmslabs.com

Robert A. Middleberg, PhD, F-ABFT, DABCC-TC, Laboratory Director

Lab ID 13000442-001

Matrix Urine

Patient Name 2492U

Patient ID 2492U

Collect Dt/Tm Not Given

Source Not Given

Analysis and Comments	Result	Units	Reporting Limit	Notes
------------------------------	---------------	--------------	------------------------	--------------

Specimens for elemental testing should be collected in certified metal-free containers. Elevated results for elemental testing may be caused by environmental contamination at the time of specimen collection and should be interpreted accordingly. It is recommended that unexpected elevated results be verified by testing another specimen in a trace metal free container.

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the US Food and Drug Administration.