



NMS Labs

CONFIDENTIAL

200 Welsh Road, Horsham, PA 19044-2208
Phone: (215) 657-4900 Fax: (215) 657-2972
e-mail: nms@nmslabs.com

Robert A. Middleberg, PhD, F-ABFT, DABCC-TC, Laboratory Director

Demo Report

Report Issued 04/27/2020 08:53
Last Report Issued 03/13/2020 12:38

88888
Clinical Example Report
Attn: Example Reports
200 Welsh Road
Horsham, PA 19044

Patient Name 8756B
Patient ID 8756B
Chain 20000432
Age Not Given DOB Not Given
Gender Not Given
Workorder 20000432
Received 03/13/2020 10:52

Sample ID 20000432-001
Matrix Blood
Patient Name 8756B
Patient ID 8756B
Container Type Clear vial

Collect Dt/Tm Not Given
Source Not Given

Approx Vol/Weight Not Given

Receipt Notes Not frozen as required

Table with 5 columns: Analysis and Comments, Result, Units, Reporting Limit, Notes

8756B Novel Psychoactive Substances (NPS)
Screen 1, Blood

Analysis by High Performance Liquid Chromatography/Time of Flight-Mass Spectrometry (LC/TOF-MS)

Table with 4 columns: Substance Name, Result, Units, Reporting Limit. Includes entries for 2-Furanylfentanyl, 3-Fluorophenmetrazine, 3-MeO-PCP, 4-ANPP, 4-MeO-PCP, 25B-NBOMe, and 25C-NBOMe.

Results for sample 20000432-001 are continued on next page



NMS Labs

CONFIDENTIAL

200 Welsh Road, Horsham, PA 19044-2208

Phone: (215) 657-4900 Fax: (215) 657-2972

e-mail: nms@nmslabs.com

Robert A. Middleberg, PhD, F-ABFT, DABCC-TC, Laboratory Director

Sample ID 20000432-001

Matrix Blood

Patient Name 8756B

Patient ID 8756B

Collect Dt/Tm Not Given

Source Not Given

Analysis and Comments	Result	Units	Reporting Limit	Notes
Synonym(s): 2C-C-NBOMe				
25H-NBOMe	None Detected	ng/mL	1.0	
Synonym(s): 2C-H-NBOMe				
25I-NBOMe	None Detected	ng/mL	1.0	
Synonym(s): 2C-I-NBOMe				
Acetyl Fentanyl	None Detected	ng/mL	0.50	
Acryl Fentanyl	None Detected	ng/mL	0.10	
Acryl fentanyl is known to have limited stability in blood which may be dependent upon pH, collection tube, and storage temperature. Negative results should be interpreted with caution.				
alpha-PVP	None Detected	ng/mL	2.0	
Synonym(s): alpha-Pyrrolidinovalerophenone; alpha-pyrrolidinopentiophenone				
Bromazepam	None Detected	ng/mL	10	
Butylone	None Detected	ng/mL	10	
Butyrylfentanyl	None Detected	ng/mL	0.10	
Synonym(s): Butyrfentanyl				
BZP	None Detected	ng/mL	10	
Synonym(s): 1-Benzylpiperazine				
Carfentanil	None Detected	ng/mL	0.10	
Synonym(s): Wildnil®				
cis-3-Methylfentanyl	None Detected	ng/mL	0.10	
Clephedrone	None Detected	ng/mL	50	
Synonym(s): 4-CMC; 4-chloromethcathinone				
Clephedrone is known to have limited stability in blood which may be dependent upon pH, collection tube, and storage temperature. Negative results should be interpreted with caution.				
Clonazolam	None Detected	ng/mL	5.0	
Cyclopropylfentanyl	None Detected	ng/mL	0.50	
Delorazepam	None Detected	ng/mL	5.0	
Synonym(s): Chlordesmethyl diazepam; Cloxazolam metabolite				

Results for sample 20000432-001 are continued on next page



NMS Labs

CONFIDENTIAL

200 Welsh Road, Horsham, PA 19044-2208
Phone: (215) 657-4900 Fax: (215) 657-2972
e-mail: nms@nmslabs.com

Robert A. Middleberg, PhD, F-ABFT, DABCC-TC, Laboratory Director

Sample ID 20000432-001

Matrix Blood

Patient Name 8756B

Patient ID 8756B

Collect Dt/Tm Not Given

Source Not Given

Analysis and Comments	Result	Units	Reporting Limit	Notes
Deschloroetizolam	None Detected	ng/mL	2.0	
Dibutylone	None Detected	ng/mL	10	
Synonym(s): bk-DMBDB				
Diclazepam	None Detected	ng/mL	20	
Ethylone	None Detected	ng/mL	10	
Etizolam	None Detected	ng/mL	10	
Flubromazepam	None Detected	ng/mL	20	
Flubromazolam	None Detected	ng/mL	5.0	
Isobutyrylfentanyl	None Detected	ng/mL	0.10	
MDPV	None Detected	ng/mL	10	
Synonym(s): Methylenedioxypropylvalerone				
Meclonazepam	None Detected	ng/mL	5.0	
Mephedrone	None Detected	ng/mL	10	
Synonym(s): 4-methylmethcathinone				
Mephedrone is known to have limited stability in blood which may be dependent upon pH, collection tube, and storage temperature. Negative results should be interpreted with caution.				
meta-Methylmethoxyacetylfentanyl	None Detected	ng/mL	0.50	
Synonym(s): 3-Methylmethoxyacetylfentanyl				
Methoxetamine	None Detected	ng/mL	2.0	
Methoxyacetylfentanyl	None Detected	ng/mL	0.50	
Methoxphenidine	None Detected	ng/mL	5.0	
Synonym(s): MXP				
Methylone	None Detected	ng/mL	10	
Methylone is known to have limited stability in blood which may be dependent upon pH, collection tube, and storage temperature. Negative results should be interpreted with caution.				
Mitragynine	None Detected	ng/mL	10	
Synonym(s): Kratom				
MPHP	None Detected	ng/mL	10	
Synonym(s): 4'-methyl-alpha-Pyrrolidinohexiophenone				

Results for sample 20000432-001 are continued on next page



NMS Labs

CONFIDENTIAL

200 Welsh Road, Horsham, PA 19044-2208
Phone: (215) 657-4900 Fax: (215) 657-2972
e-mail: nms@nmslabs.com

Robert A. Middleberg, PhD, F-ABFT, DABCC-TC, Laboratory Director

Sample ID 20000432-001

Matrix Blood

Patient Name 8756B

Patient ID 8756B

Collect Dt/Tm Not Given

Source Not Given

Analysis and Comments	Result	Units	Reporting Limit	Notes
N-Ethyl Pentylone	None Detected	ng/mL	10	
ortho-Fluorofentanyl	None Detected	ng/mL	0.10	
para-Fluorobutyrylfentanyl	None Detected	ng/mL	0.10	
Synonym(s): 4F-butyryl fentanyl				
para-Fluorofentanyl	None Detected	ng/mL	0.10	
para-Fluoroisobutyrylfentanyl	None Detected	ng/mL	0.10	
Synonym(s): 4F-Isobutyrylfentanyl; FIBF				
para-Methylmethoxyacetylfentanyl	None Detected	ng/mL	0.50	
Synonym(s): 4-Methylmethoxyacetylfentanyl				
Pentadrone	None Detected	ng/mL	2.0	
Pentadrone is known to have limited stability in blood which may be dependent upon pH, collection tube, and storage temperature. Negative results should be interpreted with caution.				
Pentylone	None Detected	ng/mL	10	
Pentylone is known to have limited stability in blood which may be dependent upon pH, collection tube, and storage temperature. Negative results should be interpreted with caution.				
Phenazepam	None Detected	ng/mL	10	
Pyrazolam	None Detected	ng/mL	5.0	
TFMPP	None Detected	ng/mL	10	
Synonym(s): 3-Trifluoromethylphenylpiperazine				
THF-F	None Detected	ng/mL	0.20	
Synonym(s): Tetrahydrofuran fentanyl				
trans-3-Methylfentanyl	None Detected	ng/mL	0.10	
U-47700	None Detected	ng/mL	1.0	
Synonym(s): U-4				
U-49900	None Detected	ng/mL	1.0	
U-51754	None Detected	ng/mL	1.0	
Valeryl Fentanyl	None Detected	ng/mL	0.50	

Results for sample 20000432-001 are continued on next page



NMS Labs

CONFIDENTIAL

200 Welsh Road, Horsham, PA 19044-2208

Phone: (215) 657-4900 Fax: (215) 657-2972

e-mail: nms@nmslabs.com

Robert A. Middleberg, PhD, F-ABFT, DABCC-TC, Laboratory Director

Sample ID 20000432-001

Matrix Blood

Patient Name 8756B

Patient ID 8756B

Collect Dt/Tm Not Given

Source Not Given

Analysis and Comments	Result	Units	Reporting Limit	Notes
------------------------------	---------------	--------------	------------------------	--------------

Sample receipt condition is inappropriate for test code: Not frozen as required

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the US Food and Drug Administration.